



12 Needles Electric Microneedle Machine 50Hz-60Hz With 6 Gear Adjustable

Our Product Introduction

For more products please visit us on hifubeauty-machine.com

Basic Information

- Place of Origin: China Guangzhou
- Brand Name: LZR
- Certification: CE
- Model Number: H5
- Minimum Order Quantity: 1
- Price: USD/60-100
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 1000 Piece/Pieces per Month



Product Specification

- Name: Electric Microneedle Machine
- Color: White
- Output Voltage: 5V
- Speed Levels: 5 Levels
- Model: M5
- Needle Number: 12 Needles
- Needle Material: Medical Grade Stainless Steel
- Input Voltage: 110-240V
- Size: 16.5cm X 4.5cm X 4.5cm
- Highlight: Electric Microneedle Machine,
6 Gear Electric Microneedle Machine,
Electric Microneedling Device 50Hz



More Images





Product Description

12 Needles Electric Microneedle Machine With 6 Gear Adjustable

Product Description:

With a frequency of 50-60Hz, this electric microneedle machine can help stimulate collagen production, reduce the appearance of fine lines and wrinkles, and improve skin elasticity. The needle number is 12, and the output voltage is 5V, making it an ideal choice for those looking for a reliable and effective microneedling device.

The Electric Micro Needle Device comes in a sleek white color, adding a touch of elegance and style to your beauty routine. Its compact design makes it easy to use at home, or even while on the go. The model number is M5, which ensures it is a high-quality machine that can deliver consistent results.

Using the Electric Micro Needle Device is simple and easy. Begin by cleaning your skin and the device, then adjust the needle length to your desired setting. Next, turn on the machine and gently glide it over your skin in a circular motion. The needles will penetrate the skin and stimulate collagen production, leaving you with a smoother, more youthful appearance.

The Electric Micro Needle Device is perfect for those looking to improve the appearance of their skin without resorting to invasive procedures. Its advanced technology and reliable performance make it an excellent choice for anyone who wants to achieve a healthier, more youthful-looking complexion.

Features:

Product Name: Electric microneedle Machine

Needle Length: 0.25mm-2.5mm

Color: White

Speed Levels: 5 Levels

Needle Number: 12 Needles

Frequency: 50-60Hz

Keywords: Cosmetic microcrystal, Beauty micro crystal pen, Beauty micro crystal pen

Applications:

The Automated Needle Therapy System is a Micro Needle Pen Device that is designed to penetrate the skin's surface, creating micro-injuries that stimulate the body's natural healing process. This device utilizes medical-grade stainless steel needles to create micro-channels in the skin, allowing for better absorption of skincare products and promoting collagen production.

With its lightweight design weighing only 1.5kg, the Automated Needle Therapy System is easy to use and maneuver, making it perfect for at-home use. It has a frequency of 50-60Hz, ensuring that it operates smoothly and efficiently.

This Micro Needle Pen Device is suitable for various skin types and can be used to achieve a wide range of skincare goals. It can be used to reduce the appearance of fine lines and wrinkles, minimize the appearance of pores, improve skin texture and tone, and even reduce the appearance of scars and stretch marks.

Whether you're looking to improve your overall skin health or target specific problem areas, the Automated Needle Therapy System is an excellent addition to your skincare routine. Use it once a week or as needed to achieve optimal results.

In summary, the Automated Needle Therapy System is a Micro Needle Pen Device that uses medical-grade stainless steel needles to create micro-channels in the skin. It has 5-speed levels, a power source of electric, weighs 1.5kg, and has a frequency of 50-60Hz. It is suitable for various skin types and can be used to achieve a wide range of skincare goals. It is perfect for at-home use and can be used once a week or as needed.

Product Parameters



Specification

Brand	Hydra pen	Model	H5
EMS intensity level	4 levels	EMS frequency	0.25 Hz
Loadable liquid capacity	3 ml	Speed Level	6 Speeds
Adjustable liquid output	3 levels	Handle Material	ABS
Needle Adjustment Range	0-1.0 mm	Adapter	US/EU/UK/AU/KR
Battery capacity	700mAh	G.W.	19.4*14.6*6.9 cm/1 pc/ 690 g
Available Catridges	H12, H36, H24, Nano-HS, Nano-HR	1 carton	41*32*38 cm/20 pcs/? kg
Customization (OEM & ODM)	logo, private label, package, new shape, change color.	Packing List	① Box ② Inductive base ③ Hydra pen H5 ④ Cable ⑤ Adapter ⑥ Manual ⑦ Retaining ring* 2 pcs ⑧ 12 pins * 2 pcs

4-pole release micro current

Using 4-pole micro current to release molecules that penetrate the bottom layer.



6-speed deep penetration

1-6 speed intelligent speed regulation solves various skin problems.



LED display screen intelligent pause

LED digital display can be used for charging and plugging, making operation more convenient.



Adjustable needle length 0~1.0 mm

Deep adjustment, The functions of regulation, repair, and hydration are fully equipped.






Serum Container:






Intuitive low battery warning

-  The number of speed flashes slowly when the battery capacity has 10% left.
-  The pen will turn off automatically after the number of speed flashes quickly for 1 minutes.
-  The pen shows "L" when restarting after the low battery automatic shutdown.


Customization:



Guangzhou Lizhirun Technology Development Co., Ltd

 18688905627

 cenglinlin25@gmail.com

 hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District, Guangzhou